

Erratum

Erratum to “Cytotoxicity evaluation of methacrylate-based resins for clinical endodontics *in vitro*”

[J Oral Sci 54, 213-217, 2012]

Esteban G. Garza¹⁾, Aniket Wadajkar²⁾, Chul Ahn³⁾, Qiang Zhu⁴⁾, Lynn A. Opperman⁵⁾, Larry L. Bellinger⁵⁾, Kytai T. Nguyen²⁾ and Takashi Komabayashi⁶⁾

¹⁾Community Dental Care Clinic, Dallas, TX, USA

²⁾Department of Bioengineering, University of Texas Arlington, Arlington, TX, USA

³⁾Department of Clinical Science, University of Texas Southwestern Medical Center at Dallas, Dallas, TX, USA

⁴⁾Division of Endodontology, University of Connecticut Health Center, Farmington, CT, USA

⁵⁾Department of Biomedical Sciences, Texas A&M Health Science Center, Baylor College of Dentistry, Dallas, TX, USA

⁶⁾Department of Endodontics, Texas A&M Health Science Center, Baylor College of Dentistry, Dallas, TX, USA

Table 1 Materials used in this study

Material name	Super-Bond C&B	Super-Bond RC Sealer	MetaSEAL	AH Plus
Abbreviation	SB-C&B	SB-RC	Meta	AH+
Manufacturer	Sun Medical/ Parkell	Sun Medical/ Parkell	Sun Medical/ Parkell	Dentsply DeTrey
Type of resin	Methacrylate resin	Methacrylate resin	Methacrylate resin	Epoxy resin
Polymerization initiator	Tributylborane (TBB)	Tributylborane (TBB)	Sulfinate salt (SS)	N/A*
Lot #	VS2	RX1	TW15 (Monomer), VG14 (Powder)	1012001095
Working time	1-2 min	6 min	30 min	4 h
Setting time	5-6 min	42 min	16 h	8 h
Clinical usage	Crown and bridge	Endodontic sealer	Endodontic sealer	Endodontic sealer

*Polymerization initiator is not necessary for polyaddition reaction.

The publisher regrets that in the above-mentioned article, Table 1 was omitted.